Colorful LED beauty mask

OPTICAL WHITENING INSTRUMENT

INSTRUCTIONS

Whitening	Freckle
Anti-age	Skin tightening

(The instruments using)

1. Clean the face of skin first. Let the water dry naturally, or use a gently washed water-absorbing towel to dry the face of the water (do not use a coarse towel to rub the face, it will lead to co ass pores or rough skin, dry skin).
2. Remove the LED color smooth mask device, then remove the power control box, connect the power supply line of the mask instrument to the output hole of the power control box.

3.according to the selection of the key to open the "Color Select / Start" the power display lights up. Displays the default time is 30 mi nut. Again by "Color Selection / Start" select the type of selection button (O 1-07 is red, blue, green, yellow, orange, purple, white), then press the "Time +", "Time -" button to set the time of "Time -" and "Time +", the simplified selection key can adjust the working time (according to individual needs, 0-60 minutes), pay attention to operation when the eyes try not to stare too much at the light source. "Brightness" can choose light intensity (O1-09).

"Flicker" can choose the light flickering mode. "Microelectron ics" can adjust micro-electric power (O1-09). (Depending on the individual needs to be adjusted, O1-09, will show the light flickering mode for the face and a comfortable and suitable position to wear, add the lightness of the mask bar, with a fixed mask. (Before wearing a mask, you can be the first to adapt to the weakest light. After wearing the finished to adjust the light intensity)

5. the photon mask can be combined with the use of disposable mask skin care products, cleaning the face after the disposable mask skin care products evenly deposited on the face, using photon irradiation, after the completion of the recommended time (S-30 minutes), the power control box automatically turns off the LED photons, at this time to remove the mask, do a good job of discoloration, so that next time

.(Working principle and function)

The high-tech LED optical whitening beauty instrument uses the principle of phototherapy. After the absorption of deep cells and better metabolism of the skin. When the light is absorbed by the skin, the light energy is converted into energy intracellular. To relax and strengthen micro blood vessels, and make the skin produce photochemical reaction-enzymatic reaction. Increase the activity of catalase (Catalase) and superoxide dismutase (SOD) in the cell, which assume the role of "guardian" and "cleaner" of the cell, increases the cellular energy source of adenosine triphosphate (ATP) decomposition. The increase in energy and protein content thus promotes the cell and the synthesis of the new replaces the old. To stimulate the secretion of basic fibroblast growth factor (BFGF) and epidermal growth factor (EGF). Collagen fibers reorganize, block and eliminate the formation of melanin. Accelerate cell growth, accelerate blood circulation, stimulate the collagen protein production, increase skin elasticity, have a great effect on repairing aging skin, acne skin, fade spots, promote and tighten, try sun burns and breastfeeding.

Efficacy of treatment

Photon treatment of 7 colors of LED lamp

• The red light (red)

• 630nm wavelength, known as "biological activity light", can enhance cell activity, accelerate cell metabolism, collagen formation at the same time, with whitening, wrinkles of skin, increase skin elasticity, improve skin condition, repair damaged skin, shrink pores, tighten skin, anti-aging antioxidant efficacy.

• Blue light (blue)

• 470nm wavelength, can inhibit inflammation, bactericidal anti-inflammatory effect, no harm to skin conditions, the protein structure is relatively damaged, in a short time to significantly reduce some inflammation to heal. Have effect on acne wound dissipation, bactericidal anti-inflammatory, repair skin, no scar and pigment.

• The green light (green)

• 520nm wavelength, with neutralization, balance and stability, to reduce skin oil secretion, balance water and oil, effectively relieve mental stress, dredge the lymph gland and swelling of the drainage.

• Yellow light (yellow)

• 580nm wavelength treatment, it can improve oxygen exchange to cells, extra energy to skin cells, decomposition of pigment, promote lymphatic drainage, improve rough skin, wrinkles, treatment of skin redness, ringworm, enhance immunity.

• Violet light (violet)

• It is the frequency of red and blue light. It combines two kinds of effective light therapy, especially has a particularly good and restorative effect in treating acne and removing acne scars. The cyanine light (cyanine)

Can gradually enhance the cell energy, the metabolism has a very good therapeutic effect.

• White light (white)

Penetrate the skin deeply, accelerate tissue metabolism active, break down spot, improve fine lines and loose skin.

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Note

1, as much as possible to avoid the light source part of the eye mask, or a long time to look at the light source. Please close your eyes and use, if you feel the light is too strong, you can block your eyes using cotton pads or a visor.

2, suggestions and disposable mask of good quality natural ingredients or the use of active whitening skin care collocation, help absorb better.

3, due to the plastic mask material, open the package is likely to have some smell, which belongs to the normal phenomenon, but the smell will soon disappear.

4, the mask is made of plastic, it is not a living light-emitting body, can not be squeezed, can not be folded, can not be weighed.

5, before using the bed, after treatment, it is best to put it back in the box, or at least get out of bed.

6, light mask can not soak in water, you can use wet cloth scrub, elastic band can be washed.

7, the mask can not be used with corrosive chemical solvent scrub.

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|Prohibit the use of the crowd|

Pregnant women, children, thyroid disease, mild allergic person, the person with skin diseases, the person with skin diseases. Skin trauma should be avoided to avoid the wound or not. Normal people use or stop using, ask the doctor doctor after doctor.

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|Product maintenance and after-sale service|

1.The company provides one year warranty service from the date of purchase, but due to time or use, the damage is included in the warranty.The product and US market parts of the company for three years, the problems occur in the product can be free repair (but product maintenance has delivery costs), more than the warranty period can also repair only change the cost of parts.

2.The fault caused by the individual, the company does not provide free warranty service: the fault caused by unauthorized disassembly, failure caused by careless beating, the resultant.

六、产品维修及售后服务(Product maintenance and after-sale service)

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七、产品维护及售后服务(Product maintenance and after-sale service)

Maintenance record(Maintenance records)

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